Homogenization Device. **PrioGENIZER™**

2006_2.00

**OPERATION MANUAL.
BETRIEBSANLEITUNG.
MANUEL D'UTILISATION.
MANUALE D'USO.
MANUAL DE INSTRUCCIONES.**

MODEL: PG003
PRODUCT NR. 10000

WWW.PRIONICS.COM

PRIONICS AG

WAGISTRASSE 27A
8952 SCHLIEREN-ZURICH
SWITZERLAND

PHONE +41 44 200 20 00
FAX +41 44 200 20 10

INFO@PRIONICS.COM
WWW.PRIONICS.COM

INDEX

1	Preface	2
2	Intended use	3
3	Safety instructions and package contents	3
4	Installation guide	4
5	Homogenization process	5
6	Description of User Interface	6
7	Operation of User Interface during homogenization process (execution mode)	6
7.1	Interrupt homogenization process	6
7.2	Stop homogenization process	6
7.3	Error and Reset	6
7.4	Device status on display	6
8	Programming revolution speed and time profiles (menu mode)	7
8.1	Revolution speed profile	7
8.2	Programming	8
8.2.1	Menu structure	9
8.2.2	Set program	9
8.2.3	Load program	11
8.2.4	Delete program	11
8.2.5	Button assignment	12
9	Cleaning and maintenance	13
9.1	PrioGENIZER™ cleaning	13
9.2	Tray cleaning	13
9.3	Rack cleaning	13
10	Repair	14
10.1	Changing a fuse	14
11	Accessories and spare parts	14
12	Error messages and trouble shooting guide	15
13	Warranty	16
14	Disposal	17
15	Appendix I: EU declaration of conformity	18
16	Appendix II: Technical Data	19
17	Appendix III: Overview	101

1] PREFACE

The Prionics PrioGENIZER™ is a homogenization device. Together with the PrioCLIP™ and the supplied racks it enables semi-automatic homogenization of samples.

You are possibly using infectious material!

For safety reasons and to ensure smooth operations you should read this manual carefully before initial operation and keep it for future reference.

GENERAL REMARKS

This operation manual contains all necessary information for the safe and correct operation and maintenance of the PrioGENIZER™. This manual is issued by the manufacturer and is part of the PrioGENIZER™. The information contained herein is intended for the operator of the device. This manual defines the intended use of the PrioGENIZER™ and contains all necessary information to safely operate it. Continuous observation of the instructions in this manual ensures the safety of people and equipment. By strictly following these instructions any manufacturer warranty will remain valid.

This manual is believed to be complete and accurate at the time of publication. In no event shall Prionics AG be liable for incidental or consequential damage in connection with or arising from the use of this manual.

Some of the information and represented particulars in photos in this manual may be different from the ones in your PrioGENIZER™. Prionics AG is continuously seeking to improve their products. The products and this manual are therefore subject to visual and technical changes without notice.

EXPLANATIONS

Important text passages in this manual are in bold type and are marked by following symbols.

NOTE

This symbol indicates situations that can result in harming human beings, material or the environment or can cause financial loss if ignored or handled incorrectly.

NOTE

This is the symbol for **biological hazard** indicating hazard to health or danger to life due to infectious material.

The term “qualified personnel” refers to staff able to conduct the necessary work when using and maintaining the PrioGENIZER™, recognizing possible hazards and avoiding them beforehand because they are experienced, well prepared and know the associated rules and regulations.

2] INTENDED USE

The main application field for PrioGENIZER™ is the homogenization of animal tissue samples. Other animal and plant material can also be processed. The use of PrioGENIZER™ for human materials is not permitted. The PrioGENIZER™ can only be used with the PrioCLIP™ designed for this purpose.

The PrioGENIZER™ processes potentially infectious materials. Work must be conducted by specially trained personnel only. National safety regulations must be strictly followed.

Do not use the device in any other way or for any other purpose because this could endanger people or equipment.

NOTE

Improper use, unauthorized modification or repair will invalidate any statutory certification of the device, or other certification which may be in force and warranty will be void.

3] SAFETY INSTRUCTIONS AND PACKAGE CONTENTS

SAFETY INSTRUCTIONS

Read the operation manual carefully before operating the PrioGENIZER™. Only trained personnel familiar with the recommended procedures for correct and safe use of the PrioGENIZER™ is allowed to operate the device.

If infectious material is processed, always wear protective clothes and gloves.

A PrioGENIZER™ contaminated by homogenate must be decontaminated and cleaned immediately (see chapter 9).

A rack contaminated by homogenate must be decontaminated and cleaned immediately (see chapter 9).

Do not operate the device if it is visibly damaged.

A defective PrioCLIP™ has to be decontaminated, sterilized and disposed of according to national safety regulations.

The PrioGENIZER™ has to be operated in a bright or well-illuminated environment in order to distinguish between processed and not yet processed PrioCLIP™.

Make sure that all PrioCLIP™ in the rack are properly closed.

Restart the homogenization process in the event of a mains power failure.

PACKAGE CONTENTS

Before installing the PrioGENIZER™ check the contents of the package to ensure no items are missing or damaged. Included in delivery:

1 PrioGENIZER™, 6 Racks, 1 Tray, 1 Power cord, 1 Manual

If the delivery is incomplete or damaged please contact your reseller.

4] INSTALLATION GUIDE

For an illustration of the PrioGENIZER™ see Appendix III.

NOTE

Homogenization can produce aerosols. Prior to work with hazardous materials a risk assessment must be conducted and national safety guidelines must be applied.

NOTE

The PrioGENIZER™ weighs approx. 25 kg and must be transported with a trolley. Two people are necessary to lift the PrioGENIZER™.

NOTE

The set-up of the PrioGENIZER™ must ensure that the main switch is accessible at all times and must allow a free space of at least 8 cm behind the ventilation slot at the back of the device, for ventilation.

Take the PrioGENIZER™ out of its packaging. Connect the PrioGENIZER™ to a grounded 100–230V AC $\pm 10\%$, 60/50Hz power outlet with the power cord provided. Do not use any adapters or multiple jacks. Jack the power cord in to the device jack (2).

Open one of the lateral flaps (6). Put the tray (7) into the device with the open side toward the display. Make sure the tray sits firmly in the device. Close the lateral flap.

Check if the limit stops (left and right, 4) protrude from the rail. If they do not, push the two triggers behind the stoppers. The limit stops do not protrude if:

- a rack was taken out of the device through the front flap
- a homogenization run was executed without a rack

Flip the main switch (1) to the **On** position (I). Close the front flap (5). The PrioGENIZER™ is now ready for use. The green **Ready** indicator light (8) is on.

5] HOMOGENIZATION PROCESS

For an illustration of the PrioCLIP™ and rack see Appendix III.

Fill the PrioCLIP™ (17) with the sample and at least 2 ml of fluid volume.
Close the PrioCLIP™ and place it in the rack (18). A maximum of 8 PrioCLIP™ can be processed at the same time.

NOTE

Make sure the PrioCLIP™ always contain enough fluid.

NOTE

Make sure that all PrioCLIP™ in the rack are properly closed.

NOTE

Make sure the limit stops (left and right) protrude. If they do not, press the two triggers behind the stoppers (see Installation guide).

1] Open the front flap (5).

2] Push the prepared rack with the filled PrioCLIP™ into the PrioGENIZER™ to the limit stops (4). The direction to push the rack into the device is marked by two arrows on the rack.

3] Close the front flap (all flaps must be closed to start the homogenization process). The green **Ready** indicator light (8) is on.

4] Press the **Start** button (13). During the homogenization process the **Ready** indicator light blinks.

5] Wait until the homogenization process is finished. The green **Ready** indicator light (8) is on.

6] Open the front flap and push the processed rack onto the tray (7).
The processed rack can also be moved to the tray by inserting the next rack.
The tray can hold up to 6 racks.

7] Close the front flap and press **Start** to process the PrioCLIP™ in the next rack.

8] If the tray is full, open a side flap (6) and take out the tray containing the processed racks. Insert an empty tray. The device cannot be started if no tray is inserted (The red **No Tray** indicator light (11) is on).

NOTE

**The last rack has to be pushed onto the tray manually.
Do not take the rack out of the device via the front flap. The limit stops do not protrude if the rack is not pushed on the tray.**

NOTE

A defective PrioCLIP™ has to be decontaminated, sterilized and disposed of according to national safety regulations.

6] DESCRIPTION OF USER INTERFACE

For an illustration of the User Interface see Appendix III.

There are two labels on every button. The upper label describes the function of the button in the execution mode [chapter 7: Operation of User Interface during homogenization process (execution mode)]. The lower label describes the function of the button in the menu mode [chapter 8: Programming revolution speed and time profiles (menu mode)].

7] OPERATION OF USER INTERFACE DURING HOMOGENIZATION PROCESS (EXECUTION MODE)

In the execution mode the upper labels of the User Interface describe the function of the buttons.

7.1] INTERRUPT HOMOGENIZATION PROCESS

Press **Stop** button (14) to interrupt the homogenization process (homogenization motors will be stopped).

Press **Start** button (13) to continue the homogenization process. The process will continue from where it was interrupted.

7.2] STOP HOMOGENIZATION PROCESS

Press **Stop** button (14) to interrupt the homogenization process (homogenization motors will be stopped).

Press **Stop** button (14) again to terminate the homogenization process (homogenization motors will be lifted).

7.3] ERROR AND RESET

The **Reset** button (15) initializes the device after an error has occurred. If the red **Error** indicator light (10) is on, an error message will occur on the upper line of the display. Look up the error message displayed in chapter 12 and execute the recommended action.

7.4] DEVICE STATUS ON DISPLAY

The current status of the PrioGENIZER™ is shown in the upper line of the display:

<i>Device status</i> <i>Name of current program</i>
--

[FIGURE 1] Display in execution mode.

[TABLE 1] Overview of device status.

Device status	Indicators				Remarks
	Ready	Stop	Error	No tray	
Initializing...	Off	Off	Off	Off	Device is initializing (after power on).
Ready	On	Off	Off	Off	Device is ready. Menu mode can be activated or homogenization process can be started.
Coupling...	Blinking	Off	Off	Off	The homogenization motors and the PrioCLIP™ are coupled.
Homogenizing...	Blinking	Off	Off	Off	Samples are being homogenized.
Decoupling...	Blinking	Off	Off	Off	The homogenization motors and the PrioCLIP™ are decoupled.
Stopped	Blinking	On	Off	Off	A started process was interrupted by pressing the Stop button. The Ready indicator blinks to show the process was interrupted.
Left/Front/Right flap open/ flaps open	On/Off/ Blinking	On	Off	Off	One or more flaps were opened. The Ready indicator stays in its previous state.
No Tray	On/Off/ Blinking	On	Off	On	No tray in the device. The Ready indicator stays in its previous state.
Error	Off	Off	On	Off	An error has occurred. Press the Reset button to confirm.
Initializing menu...	On	Off	Off	Off	Menu mode is initializing.
Menu	Off	Off	Off	Off	Device is in menu mode.

8] PROGRAMMING REVOLUTION SPEED AND TIME PROFILES (MENU MODE)

The characteristics of the revolution speed significantly influence the homogenization result. 9 different speed-time characteristics can be defined with the PrioGENIZER™. The factory set program **P0 PRIONICS TSE** cannot be changed or deleted.

8.1] REVOLUTION SPEED PROFILE


The profile of the speed is defined through a curve containing a maximum of 12 points (see Figure 2: Speed profile).

Each point of the curve consists of the two parameters *speed and time*.

[TABLE 2] Range of parameters to define a speed profile.

Range for <i>speed [rpm]</i> :	0....24'000 rpm in steps of 500 rpm
Range for <i>time [t]</i> :	0....99s in steps of 1s

The revolution speed is set absolutely while the time is set relatively for each parameter. The time of one homogenization cycle is the sum of all the time parameters. The speed for parameter 12 is always 0.

[FIGURE 2] Speed profile.

8.2] PROGRAMMING

To enter into the programming mode press **Menu** (16) while the device is in the ready status. In the menu mode the lower labels indicate the function of the buttons. The main menu “Select menu” appears on the display after entering the menu mode.

Select menu
Set program

[FIGURE 3] Display after entering the menu mode.


Navigate with the **Down** (14) or the **Up** button (13) through the main menu. The following sub menus can be chosen:

[TABLE 3] The sub menus “Set program”, “Load program” and “Delete program” can be selected in the main menu. Select “Quit program” to exit the menu mode.

Set program	Program new revolution speed profile.	Chapter 8.2.1
Load program	Load revolution speed profile already programmed.	Chapter 8.2.3
Delete program	Delete revolution speed profile already programmed.	Chapter 8.2.4
Quit menu	Quit menu.	

- Press **Confirm** (15) to select the indicated menu.
- Press **ESC** (16) to quit the menu mode.

8.2.1] MENU STRUCTURE

[FIGURE 4] **Menu structure.**

8.2.2] SET PROGRAM

Select “Set program” from the main menu. 9 different programs can be set (P1 to P9). If a program has already been defined, the selected program name appears behind the program number. If no parameters are defined under the selected program number, the display reads “<< EMPTY >>”.

```
Set program
P1 <<EMPTY>>
```

[FIGURE 5] **Program number 1 is chosen to set new parameters.**

- Select the program number by pressing **Up** or **Down**.
- Confirm the program number by pressing **Confirm**.

First enter the program name

Scroll through the alphabet using **Up** and **Down** (Character sequence: [“Space”, A, B, C, ... Z, 0, 1, ...9]). Select the displayed letter by pressing **Confirm**. By pressing **Confirm** the next letter can be selected.

```
Set program name
P1 TES
```

[FIGURE 6] **Display when entering a program name.**

- Press **Up** to select the next alphabetic character or hold **Up** for more than 1 second to scroll forward through the alphabet more rapidly.
- Press **Down** to select the previous alphabetic character or hold **Down** for more than 1 second to scroll backward through the alphabet more rapidly.
- After completing the program name hold **Confirm** for at least 2 seconds to confirm the program name.

Enter the time parameter

```
Set program param1
t: 0 rpm: 0
```

[FIGURE 7] Display when entering the first time parameter.

- Press **Up**, the time value increments by 1 s. Hold **Up** for more than 1 second to increase the time value in steps of 10 s.
- Press **Down** to decrease the time value by 1 s. Hold **Down** for more than 1 second to decrease the time value in steps of 10 s.
- Press **Confirm** when the desired time value is displayed.

The maximum time value is 99 s.

Enter the speed parameter

```
Set program param1
t: 15 rpm: 0
```

[FIGURE 8] Display when entering the first speed parameter.

- Press **Up** to increase the speed by 500 rpm. Hold **Up** for more than 1 second to increase the speed in steps of 1000 rpm.
- Press **Down** to decrease the speed by 500 rpm. Hold **Down** for more than 1 second to decrease the speed in steps of 1000 rpm.
- Press **Confirm** when the desired speed value is displayed.

The maximum speed value possible is 24'000 rpm.

```
Set program param2
t: 0 rpm: 0
```

[FIGURE 9] Display when entering the second time parameter.

Repeat the steps as described above until the entire speed profile is defined. A maximum of 12 parameters can be defined. The speed of the 12th parameter is always 0 and cannot be changed.

After entering the last parameter hold **Confirm** for at least 1 second. The remaining parameters are set to zero and the following message is displayed:

```
Save program [Co/Es]
"name"
```


[FIGURE 10] Message to save a program. On the lower line the program name is displayed.

- Press **Confirm** to save the program.
- Press **ESC** to return to the main menu without saving.

The "Set program mode" can be quit at any time. Hold **ESC** for at least 1 second to go back to the main menu without saving the program.

8.2.3] LOAD PROGRAM


Select “Load program” from the main menu. Navigate with the **Down** or the **Up** button through the main menu. Press **Confirm** when “Load program” appears on the lower display line.

```
Select menu
Load program
```

[FIGURE 11] Display in the main menu mode to select the “Load program” sub menu.

10 different programs can be loaded (P0 to P9). P0 is the factory set program “PRIONICS TSE”. P1 to P9 can be defined by the user. When a program is defined, the program name appears behind the program number. If no parameters are defined under the selected program number the display reads “<< EMPTY >>”.

```
Load program
P0 PRIONICS TSE
```


[FIGURE 12] Factory set program number 0 is chosen to load.

- Navigate with the **Down** or the **Up** button through the different program numbers.
- Press **Confirm** to load the selected program.
- Press **ESC** to exit the sub menu.

The next homogenization process can be performed with the selected program.

8.2.4] DELETE PROGRAM


Select “Delete program” from the main menu. Navigate with the **Down** or the **Up** button through the main menu. Press **Confirm** when “Delete program” appears on the lower display line.

```
Select menu
Delete program
```

[FIGURE 13] Display in the main menu mode to select the “Delete program” sub menu.

The factory set program “P0 PRIONICS TSE” cannot be deleted.

```
Delete program
P1 TEST1
```

[FIGURE 14] Program number 1, called TEST1, is chosen to be deleted.

- Navigate with the **Down** or the **Up** button through the different program numbers.
- Press **Confirm** to delete the selected program.
- Press **ESC** to exit the sub menu.

The program name and all parameters are deleted. The deleted program number will appear as <<EMPTY>>. The next homogenization process will be performed with the currently selected program. If the loaded program was deleted, “P0 PRIONICS TSE” is selected by default.

8.2.5] BUTTON ASSIGNMENT

[TABLE 4] Button assignment in “Set program”.

Button	Action when pressed	Action when held >1 s
Start/Up	Backwards	–
Stop/Down	Forwards	–
Reset/Confirm	Select sub-menu	–
Menu/ESC	Quit sub-menu	–

[TABLE 5] Button assignment in “Load program”.

Button	Action when pressed	Action when held >1 s
Start/Up	Backwards	–
Stop/Down	Forwards	–
Reset/Confirm	Load program	–
Menu/ESC	Quit sub-menu	–

[TABLE 6] Button assignment in “Delete program”.

Button	Action when pressed	Action when held >1 s
Start/Up	Backwards	–
Stop/Down	Forwards	–
Reset/Confirm	Delete program	–
Menu/ESC	Quit sub-menu	–

[TABLE 7] Button assignment in “Set program name”.

Button	Action when pressed	Action when held >1 s
Start/Up	Next character	Fast forward
Stop/Down	Previous character	Fast backward
Reset/Confirm	■ Next position ■ At the “ok” position: Change to menu “Set program param1”	Save and return to menu “Set program”
Menu/ESC	■ Previous position	Return to menu “Set program” without saving

[TABLE 8] Button assignment in “Set program param1”.

Button	Action when pressed	Action when held >1 s
Start/Up	■ Increases time value by 1 s ■ Increases speed value by 500 rpm	■ Increases time value by 10 s ■ Increases speed value by 1000 rpm
Stop/Down	■ Reduces time value by 1 s ■ Reduces speed value by 5000 rpm	■ Reduces time value by 10 s ■ Reduces speed value by 1000 rpm
Reset/Confirm	■ Next position ■ At the “ok” position: Change to next parameter	Save and return to menu “Set program”
Menu/ESC	■ Previous position ■ Back to previous parameter	Return to menu “Set program” without saving

9] CLEANING AND MAINTENANCE

The PrioGENIZER™ processes potentially infectious materials. Cleaning must be conducted by specially trained personnel only. National safety regulations must be strictly followed.

A PrioGENIZER™ contaminated with homogenate must be decontaminated and cleaned immediately.

A rack contaminated with homogenate must be decontaminated and cleaned immediately.

Do not operate the device if it is visibly damaged.

A defective PrioCLIP™ has to be decontaminated, sterilized and disposed of according to national safety regulations.

9.1] PrioGENIZER™ CLEANING

Clean the PrioGENIZER™ with 70% ethanol solution.

NOTE

Always wear protective clothes and gloves.

NOTE

Always unplug the mains before cleaning.

NOTE

Many cleaning agents are easily inflammable. Do not smoke while using cleaning agents. Open flames are prohibited.

9.2] TRAY CLEANING

To decontaminate the tray (7) use 2M sodium hydroxide (2M NaOH).

To clean the tray use 70% ethanol solution.

9.3] RACK CLEANING

To decontaminate the rack (18) use 2M sodium hydroxide (2M NaOH).

To clean the rack use 70% ethanol solution.

10] REPAIR

For an illustration of the position of the fuse compartment see Appendix III.

10.1] CHANGING A FUSE

To change a fuse in the PrioGENIZER™, proceed as follows:

- 1] Switch off the device with the main switch (1) (Position "0").
- 2] **Disconnect the power cord from the mains outlet and unplug it from the PrioGENIZER™ (2).**
- 3] Open the fuse compartment (3) below the jack using a screwdriver.
- 4] Put a replacement fuse type FST5 x 20/4A/250V Time-lag in the designated compartment (3).
- 5] Close the compartment.
- 6] Jack the power cord into the mains outlet and the device jack (2).
- 7] Switch-on the device by pressing the main switch (1) (position "I"). If all flaps are closed, the tray is properly inserted and the initialization of the device was successful, the green **Ready** light (8) will be on.

NOTE

For any other repairs contact your reseller. Prionics AG cannot be made responsible for damages caused by incorrect or unauthorized repairs.

11] ACCESSORIES AND SPARE PARTS

PrioCLIP™ Order No. 10010.

Container to hold tissue sample and liquid solution.

Rack Order No. 10070.

Holder for 8 PrioCLIP™.

Tray Order No. 10060.

Receptacle for maximum 6 racks.

Storage System Order No. 10080.

Trolley for storage and transport of up to 6 trays.

12] ERROR MESSAGES AND TROUBLE SHOOTING GUIDE

Following errors can occur during operation of the device. If an error occurs, the red **Error** light is on. The error message appears on the upper line of the display.

[TABLE 9] **Error messages.**

Problem	Can occur in status	Possible cause	Action
Error: Overcurrent CM	Initializing... Coupling... Decoupling...	An over-current in the gear motor was detected. Cause: ■ Rack got jammed ■ Rack was not sitting firmly on the limit stops ■ Coupling of the PrioCLIP™ failed (stirring rod possibly damaged)	Press Reset . If the error message recurs contact your reseller.
Error: Timeout CM	Initializing... Coupling... Decoupling...	The beam with the homogenization motors did not reach the coupled or de-coupled position within the designated time. Cause: ■ Gear motor defect ■ Sensor defect	Press Reset . The device initializes. If the error recurs contact your reseller.
Error: MB x Motor y	Coupling... Homogenizing...	One of the 8 homogenization motors did not reach its designated speed. MB x (x=1 or 2) indicates on which motor board the error occurred and Motor y (y=1..4) indicates which motor did not reach the speed. Error cause: ■ Coupling of the PrioCLIP™ failed (stirring rod possibly damaged) ■ Homogenization motor possibly out of order	Press Reset and start homogenization process again pressing Start . The device initializes. If the error message recurs contact your reseller.
Error: SW HP, Error: HP, Error: SW CP, Error: CP, Error: SW SM	Coupling... Homogenizing... Decoupling...	These messages indicate unusual software-errors.	Press Reset to initialize the device. Report the error message to your reseller.

[TABLE 10] Trouble shooting guide.

Problem	Can occur in status	Possible cause	Action
Homogenization does not start.	Green Ready light is on.	Start was not pressed.	Press Start .
Homogenization does not start.	Red Stop light is on.	The device is not ready.	Check the status readout in the first line of the display (e.g. close the flaps, insert tray).
Homogenization stops, Ready indicator blinks	Red Stop light is on.	A running process was interrupted by pressing Stop .	Restart the process by pressing Start again.
Homogenization stops	Red Error light is on.	An error occurred during a homogenization run.	See Table 9. Error messages.
	No light is on.	Various causes.	Execute following steps: ■ Check the power cord is properly connected ■ Check the main switch is on ■ Check the fuse ■ Call your reseller

13] WARRANTY

Damages caused by abnormal use, transport, unauthorized modifications or misuse, improper or inadequate maintenance of the device are excluded from the warranty. Damages caused by material faults or production errors will be repaired or the defective part replaced by the manufacturer free of charge. The Prionics AG commercial warranty is 1 year.

Customers may have other additional legal rights which vary from country to country and which are not affected in any way by this commercial warranty.

Prionics AG will not provide warranty repair or replacement if in our opinion the problem resulted from externally caused damage or use outside of the product specification or from the use of options, parts or consumables which are not Prionics branded or approved.

In no event will Prionics AG be liable for special direct or indirect or consequential damages, including loss of profits or inability to use the product, even if Prionics AG has been advised of the possibility of such damages. In no event shall Prionics AG liability for any damages or loss to the current user or any other party exceed the price paid for the product.

14] DISPOSAL

The legal disposal of the PrioGENIZER™ is the responsibility of the purchaser. You may contact your dealer for further details.

NOTE

A defective PrioCLIP™ has to be decontaminated, sterilized and disposed of according to national safety regulations.

NOTE

During its use the PrioGENIZER™ was possibly contaminated with potentially infectious materials. Please follow the following instructions carefully!

You must follow this procedure before the PrioGENIZER™ can be disposed of professionally:

- 1]** Switch off the PrioGENIZER™
 - 2]** Unplug the power cord
 - 3]** Decontaminate the PrioGENIZER™ according to local national safety regulations.
-

NOTE

In case of difficulties concerning the disposal you may contact your reseller!

15] APPENDIX I: EU DECLARATION OF CONFORMITY**ENGLISH**

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
16] APPENDIX II: TECHNICAL DATA

Model name	PrioGENIZER™
Operation voltage	100–230V AC $\pm 10\%$, 60/50Hz
Power consumption	85 W (max)
Fuse	5x20, 4A/250V/T
Max. relative humidity	80% at 30°C
Max. operation altitude	2000 AMSL
Temperature	0°C to 45°C
Noise pressure level [L _{pA,1m}]	73 dB(A)
Dimensions [mm]	450 x 500 x 300
Weight	approx. 25 kg
SW-Version	Press Reset/Confirm in Ready mode and the SW-Version will be displayed briefly
Pollution Degree	2
Overvoltage category	II

Transport and storage condition

Storage temperature	0°C to 70°C
Transport temperature	–10°C to 70°C

- The identification and LOT-No. is placed at the back side of PrioGENIZER™.
- Subject to technical changes and error.

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17] APPENDIX III: OVERVIEW

17] APPENDIX III: OVERVIEW

1	MAIN SWITCH On / Off
2	JACK For power cord.
3	FUSE BOX
4	LIMIT STOPS Left and right in rack slot.
5	FRONT FLAP
6	FLAPS LEFT / RIGHT
7	TRAY
8	READY INDICATOR Ready for operation / Operating (blinking green).
9	STOP INDICATOR Homogenization stopped (permanent red).
10	ERROR INDICATOR Error occurred (permanent red).
11	NO TRAY INDICATOR Tray not present (permanent red).
12	DISPLAY Shows operating status and program parameters.
13	START/UP BUTTON Starts Homogenization process. Used for programming.
14	STOP/DOWN BUTTON Interrupts Homogenization process. Used for programming.
15	RESET/CONFIRM BUTTON Resets PrioGENIZER™ process to original status. Confirms Menu entries.
16	MENU/ESC BUTTON Starts programming menu. Cancels menu entries.
17	PrioCLIP™
18	RACK

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